



A General Protocol for PEIone™ Transfection Reagent

9620 Medical Center Dr, Ste 100, Rockville, MD 20850

Tel. +1 301-605-7502

Email: customer@aavnergene.com

Website: www.aavnergene.com

This Product is for Research Use Only -

Not for Diagnostic or Treatment Purposes

SKU: AR000000

Product Description

PEIone™ is a linear polyethylenimine (PEI)-based transfection reagent specifically optimized for high-efficiency adeno-associated virus (AAV) production in suspension and adherent HEK 293T or HEK 293 cell systems. As a chemically synthesized, animal-origin-free reagent, PEIone™ ensures low cytotoxicity while maintaining robust and reproducible transfection performance.

Supplied as a clear, sterile, ready-to-use liquid, PEIone™ is fully compatible with our novel AAVone® single-plasmid system, and other plasmid-based AAV production systems, such as our AAVdual® two-plasmid system and traditional AAVtri three-plasmid system, making it broadly applicable in both research and large-scale AAV manufacturing.

Storage

- Store at 2–8 °C
- Stable for up to 2 years
- Do not freeze

Key Advantages

- High DNA-binding capacity with efficient performance at low reagent-to-DNA ratios
- Uses only one-third of the reagent, reducing consumption and overall cost
- Stable complex formation ensuring consistent transfection efficiency
- Scalable for both research- and industrial-scale AAV production

Protocol

1. Required Materials

- HEK 293 or HEK 293T cells (suspension or adherent).
- Plasmid DNA (pAAVone® plasmids or other plasmids for AAV production; 0.5–5 µg/µL stock).
- Appropriate cell culture medium.

2. Preparing Cells

- Suspension cells: Grow cells to a density of $4.0\text{--}6.0 \times 10^6$ viable cells/mL with $\geq 95\%$ viability. Dilute to $2.5\text{--}3.5 \times 10^6$ viable cells/mL using pre-warmed medium at the time of transfection.
- Adherent cells: Seed cells to reach 70–90% confluence at the time of transfection.

3. Transfection Protocol for 240 ml culture

- Bring PEIone™ and culture medium to room temperature.
- Dilute plasmid DNA in culture medium to ~5% of the culture volume (e.g., 360 µg DNA in 12 mL medium for 240 mL culture). Mix gently.
- In a separate tube, dilute PEIone™ reagent in culture medium to ~5% of the culture volume (e.g., 126 µL PEIone™ reagent in 12 mL medium for 240 mL culture). Mix gently.
- Add diluted PEIone™ solution to the diluted DNA solution. Mix gently.
- Incubate for 15–30 minutes at room temperature to allow complex formation.
- Add complexes directly to the culture vessel.
- Incubate cells for about 68-72 hours, then harvest according to in-house AAV production protocols.

Recommended DNA and Reagent Amounts

The amount of PEIone™ and DNA, and their ratios should be optimized for your AAV genomes and serotypes/capsids to obtain maximum productivity.

Typical PEIone™/ pAAVone® Ratios:

- Suspension HEK 293 or HEK 293T cells: PEIone™/DNA = 0.5:1–1:1
- Adherent HEK 293 or HEK 293T cells: PEIone™/DNA = 1:1–2:1

Typical pAAVone® amount:

- Suspension HEK 293 or HEK 293T cells: 0.3-0.6 µg per 1×10^6 cells
- Adherent HEK 293 or HEK 293T cells: 0.3-0.6 µg per well (24-well plate)

Typical PEIone™ amount:

- Suspension HEK 293 or HEK 293T cells: 0.2-0.5 µl per 1×10^6 cells
- Adherent HEK 293 or HEK 293T cells: 0.5-1.0 µl per well (24-well plate)

Table 1. Recommended DNA amounts and reagent volumes for transfection of pAAVone® plasmid at various scales in suspension HEK 293/293T cells

	Vessel type				
	6-well	125 mL flask	250 mL flask	500 mL flask	1 L flask
Cell number	6×10^6	90×10^6	180×10^6	360×10^6	720×10^6
Culture volume	2 mL	30 mL	60 mL	120 mL	240 mL
DNA amount	3 µg	45 µg	90 µg	180 µg	360 µg
PEIone™	2.1 µL	31.5 µL	63 µL	126 µL	252 µL
Dilution medium	2×100 µL	2×1.5 mL	2×3 mL	2×6 mL	2×12 mL

Table 2. Recommended DNA amounts and reagent volumes for transfection of pAAVone® plasmid at various scales in adherent HEK 293/293T cells.

	Vessel type					
	24-well	12-well	6-well	T-25 flask	T-75 flask	100 mm dish
Culture medium	0.5 mL	1 mL	2 mL	4 mL	10 mL	10 mL
DNA amount	0.5 µg	1 µg	2 µg	4 µg	8 µg	8 µg
PEIone™	0.75 µL	1.5 µL	3 µL	6 µL	12 µL	12 µL
Dilution medium	2×25 µL	2×50 µL	2×100 µL	2×200 µL	2×0.5 mL	2×0.5 mL

Notes

For triple-plasmid AAV systems (AAVtri), the total plasmid amount required is typically higher than that for pAAVone®. The ratio of the three plasmids should be optimized for maximum AAV productivity. A starting molar ratio of 1:1:1 for each plasmid is recommended for initial testing.